



***Optimizing Sampling and Analytical Parameters
for Soil Vapour Samples using Automated Thermal
Desorption / Gas Chromatography / Mass
Spectrometry (ATD/GC/MS)***



Stephen Varisco, Technical Manager, CARO Analytical Services
Lee Marotta, GC/MS Product Specialist, PerkinElmer LAS
Miles Snow, GC/MS Product Specialist, PerkinElmer LAS

Outline

Objective: Soil Vapour Tube

➤ Background

- **CARO Analytical Services**
- **Soil Vapour Regulations**

➤ Thermal Desorption Technology

- **Advantages**
- **Operation**

➤ Tube Design and Adsorbent Optimization

- **Limitations of current tubes for soil vapour samples**
- **Objectives of new tube**

➤ New Tube Performance Characteristics

Background information on CARO Analytical Services

- **British Columbia based**
 - **Full Service Environmental Testing Capabilities**

- **Technical Leadership in Soil Vapour**
 - **Co-authored BC SV Reference Methods**
 - **ISO 17025 Accreditation for SV (www.cala.ca)**

- **Client Collaboration**
 - **Knowledgeable & Accessible Staff**
 - **Equipment Rentals, Client Training**

- **SR&D Projects**
 - **VOC Artifacts in Tubing (www.caro.ca/soilvapour)**
 - **Custom TD Tube Development**
 - **VOC Sorption on Tubing – in progress**

Regulatory Background

- **23 US States regulate Soil Vapour Intrusion**

- **British Columbia Contaminated Sites Regulations, Schedule 11**
 - **Effective January 1, 2009**
 - **118 regulated parameters: >90% via TD**
 - **Stringent regulations**
 - **Trichloroethene (TCE) – 1.0 ug/m³ (dry cleaning site)**
 - **Benzene – 1.5 ug/m³ (gasoline site)**

- **Why is soil vapour important?**
 - **Protection of Human Health**

- **Test water and soil first? Or ... Test vapour first?**



Advantages and Operation of Thermal Desorption Technology


Advantages of Thermal Desorption Technology

Air Toxics ≠ Charcoal Tubes


- 1. Established methodology**
- 2. Convenient transport**
- 3. Easy to clean ... immediate reuse ... fast turnaround**
- 4. Cost effective**
- 5. Large sample volumes**
- 6. Suitable for non-polar and polar compounds**
- 7. Inherent Water Management**
- 8. Enables Recollection to preserve sample**



State 1: Sample Tube Desorption



Stage 2: Transfer of Sample to Instrument



Optimizing Sample Integrity with Thermal Desorption

- **Automatic Internal Standard addition**
- **Recollection onto the same tube or new tube**
- **Automated Leak Checking**
- **Impedance Testing**
- **Excellent Water Management**




Total Ion Chromatogram of a 74 TO-15 and TO-17 Component Mix

Linearity

AT_032609_09

27-Mar-2009 + 01:37:40

Scan E+
TIC
3.23e8






Tube Design and Adsorbent Optimization

Thermal Desorption Tube

Multiple Adsorbents

- accommodate wide boiling point range sample



Goal – Optimize a Tube for Soil Vapour Intrusion Sampling

- **Soil vapour differs from other air sampling applications:**
 - **Higher moisture**
 - **Greater analyte range**
 - **Wider concentration range**

- **Limitations of current industry-standard sorbent tubes:**
 - **Lack of a weak adsorbent prevents release beyond nC12**
 - **Ensure Safe sampling volumes for the low boilers**

Characteristics of New Thermal Desorption Tube

- **Broadest Analytical Range**
 - dichlorodifluoromethane to phenanthrene
 - nC3 – nC22+

- **Protects the Strong Adsorbents**
 - Prevents irreversible adsorption
 - Clean after one desorption cycle

- **Good Safe Sampling Volumes**

- **Optimal Water Management**



Soil Vapour Intrusion Tubes Performance Characteristics

Precision, Linearity (average for Class) and Reporting Limits



Instrumentation:

PerkinElmer TurboMatrix 650 Thermal Desorber


PerkinElmer Clarus 600 GC/MS

Compound Class	# Comp	Precision (%RSD)	Dynamic Range	Reporting Limit	Outliers
		n=8	0.2 to 200 ng	10 liter sample (FS)	on RL
Gasses	6	6.9%	0.9952	0.05 ug/m3	
non-Aromatic Halogens	33	2.7%	0.9985	0.02 ug/m3	CHCl ₃ =0.05ppb
Aromatics	15	1.4%	0.9995	0.02 ug/m3	
Halogenated Aromatics	9	1.4%	0.9997	0.02 ug/m3	
Others	10	2.7%	0.9965	0.05 ug/m3	

Recovery Results EXCELLENT

- **79 components investigated**
 - 400ng VOC mix + 250ng PAH
 - Range: difluorodichloromethane to phenanthrene
- **Recovery procedure**
 - Analyzed spiked tube
 - Analyzed blank tube
 - Re-analyzed spiked tube which should be clean
- **Non-detectable carryover**
 - Slight carryover of 4 heaviest PAHs
 - Significantly below method criterion

PAH Compounds	% Carryover
	(very low)
1-Methyl Naphthalene	0.3
Anthracene	0.2
Fluorene	0.6
Phenanthrene	1.2



Tube Robustness – Diesel


Carryover <1%

08-Apr-2009 + 16:57:46

new 650_040809_08

Scan EI+
TIC
6.96e8


10,000ng Diesel spiked tube



new 650_040809_10

Scan EI+
TIC
6.96e8


Reanalysis of this tube to determine Carryover



What is Breakthrough?

- **Goal: Regulated compounds do not break through**
- **EPA TO-17:** “The volume sampled when the amount of analyte collected in a back-up sorbent tube reaches a certain percentage (typically 5%) of the total amount collected by both sorbent tubes”
- **Calculation: $BT / (BT + \text{spiked}) * 100\%$**
 - **BT = concentration of back tube**
 - **Spiked = concentration of front tube**
- **Safe Sampling Volume (SSV): 2/3 of Breakthrough Volume**
- **Assessing breakthrough in the real-world:**
 - **Sampling in series**
 - **Sampling in parallel (distributed volume pairs)**


Manifold used for determining breakthrough & moisture retention



Sample Tube spiked

- 300ng of 79 components
- Total = 23,700ng on tube
- Passed humidified N₂ through tube
 - flow rate = 100mL/min
 - Time = 100min
 - Mimic 10L sampling
 - All studies performed in triplicate


Slight Breakthrough of only two Components out of 79 VOCs



VOC gasses:

- 10L of sample
- 70% humidity

Component	% BT
Dichlorodifluoromethane	1.0
Chloromethane	5.4
Vinyl Chloride	nd
Bromomethane	nd
Chloroethane	nd
Trichlorofluoromethane	nd



Water Management

Methods of Water Removal

1. Nafion Drier / Desiccants

- **Polar Compounds Removed - Cannot be used**

2. Hydrophobic adsorbents


3. Minimize sampling volumes while maintaining regulated detection limits

4. Dry Purging

- **Typically not required for indoor air**


Why Remove Moisture?

- Mass Spectrometer (detector) problems
 - Signal quenching
 - Increased maintenance
- Poor chromatography
 - Can cause false negatives



Real World Sample – Petroleum Contamination (Tank Pull Site)

Component	Result (ug/m3)
1,2,4-trimethylbenzene	72.7
1,3,5-trimethylbenzene	16.2
Benzene	1.6
Dichlorodifluoromethane	9.1
Ethylbenzene	12.3
Naphthalene	20.9
Tetrachloroethene	9.9
Toluene	45.7
Trichlorofluoromethane	1.2
Xylenes (total)	80.0
nC6-nC13	9950



Summary

- **Thermal Desorption Technology**
 - **Tubes well-suited for vapour intrusion investigations**
 - **Instrumentation Advancements → Analytical Integrity**

- **Team developed new Thermal Desorption Tube that Achieves**
 - 1. Broadest Analytical Range**
 - 2. Protection of the Strong Adsorbents**
 - 3. Good Safe Sampling Volumes**
 - 4. Optimal Water Management**

- **Presentation available for download at www.caro.ca or www.perkinelmer.com**

Acknowledgements

Project Team:

➤ CARO Analytical Services:

- **Stephen Varisco, Technical Manager**
- **Luba Tsurikova, Senior Analyst**
- **Patrick Novak, Business Manager**
- **Brent Mussato, President**

➤ PerkinElmer:

- **Lee Marotta, Senior Product Specialist**
- **Miles Snow, Senior Product Specialist**